

UMS DYNAMIC UFM



Urodynamic modular system DYNAMIC UFM

is intended for a urodynamic examination of lower urinary tract. A record evaluation enables to the doctor to determine a dysfunction of urethra or detrusor, to choice the adequate treatment method according to the incontinence type. Urodynamic system is suitable both for gynecological and urological ambulances or hospital facilities.

For data displaying and processing, system DYNAMIC utilizes a classical computer with operating system MS Windows.

DYNAMIC UFM enables measuring and displaying of measured quantities, their processing and evaluating by the form of graphs and tables and enumerations of parameters. The relevant results are printed in color by the form of graphs and tables and archived on the computer together with personal data of patient.

up to **7**
years warranty

UMS DYNAMIC UFM



- + Advanced ergonomic design
- + Stainless steel UFM stand
- + Comunication by cable or wireless